

1 Person Articulating Decon Shelter



Purpose

An individual decontamination shelter for effective decontamination of a single ambulatory person.

Features & Benefits

- For effective decontamination of single ambulatory person.
- Shelter incorporates a tall profile to accommodate individuals fully outfitted in PPE
- Patented external articulating frame with corrosion resistant anodized aluminium
- Design facilitates rapid deployment
- Shelter can operate with 30% of frame missing or broken
- UV stable, fire retardant and chemical resistant fabric on interior. Provides durability; protects exterior frame from contaminants and expedites post-use clean up
- Space efficient for facilities with limited space
- Interior Bund covers entire floor of decon area; contains the waste water run-off
- On/Off capable shower nozzles on interior fabric saturate patients with 360° of spray



Technical Data

- Deployed Dimensions(m): L 1.50 x W 2.1 x H 2.55
- Internal Dimensions (m): L 1.5 x W 1.7 x H 2.3
- Packed Dimensions(m): L 0.36 x 0.86 x 0.81
- Weight: 43kg
- Shelter Set Up: 1 people, 2 minutes
- Shelter Strike Down: 1 people, 2 minutes

Quantity	Description
1	Decontamination shelter (with integrated shower system)
1	Interior berm to contain water run-off
7	Integrated shower nozzles
2	Zippered panel end doors (1 each end)
1	Hose adaptor kit
1	Anchor kit
1	Shelter Repair kit
1	Transport Bag



Unit 29 Arkwright Road, Astmoor Industrial Estate, Runcorn, Cheshire, WA7 1NU
 Tel: +44(0)1928 575307 Email: info@grs-uk.com Web: www.grs-uk.com